

3M is boldly redefining
the economics of
CAD/CAM dentistry ...



putting it
within reach
for every dentist.

3M

True Definition Scanner

Predictable clinical outcomes

Benefits of digital process versus traditional process

Our world engages with the speed and unlimited capabilities of digital. For the dental practice, that means not only faster, easier patient care, but better care. Digital impressions eliminate inaccuracies that are inherent to the traditional process, making clinical outcomes more predictable.

Better patient experience

Digital technology transforms the procedure into an interactive educational experience. Seeing images on screen gives patients a true understanding of their oral situation and what is needed.

It's easy to see why dentistry is going digital.


Better oral care

With powerful tools that are unique to digital dentistry, you can view each case in amazing detail and assess your work on the spot, eliminating many inaccuracies. Plus, digital files enable better communication with your lab for better results.

Improved productivity

Digital processes elevate dentistry for you and your staff in ways that can have a profound impact on your practice.

- Less stress for you and your staff
- Fewer procedural steps and less hassle
- Faster turnaround and more accurate results from the lab
- Reduced seating times and remakes



Digital Impressions (Model Free) Workflow. A model is no longer necessary. The accuracy of the data ensures that the designer is operating on a precise replica of the oral anatomy. Model-free digital impressions deliver an accurate, quality restoration while saving on time and lab labor.



The 3M™ True Definition Scanner provides multiple workflow options. Learn more in the "How it Works" section at www.3M.com/TrueDefinition



True Definition Scanner

More accurate—and more consistently accurate—than leading systems on the market.¹

The 3M™ True Definition Scanner features the powerful “3D-in-motion” video technology, which allows you to capture and simultaneously view a true replica of the oral anatomy in extraordinary detail. This gives the flexibility needed for quadrant dentistry and the accuracy needed for full arch dentistry. Now, this advanced technology is enhanced by a small and comfortable wand and interactive multi-touch monitor.

Choices supported by an open and secure cloud-based platform to store, share and connect

The 3M™ Connection Center offers more choices for sharing clinical case files—including open STL files for flexibility. You, your dental lab or other partners have the flexibility to choose the best solution for the desired clinical outcome.

The 3M True Definition system is designed to connect with mills and other systems—current and future—making it a smart investment for any practice.



Small and robust. Designed without fragile moving parts, so it never needs calibrating. The scanner is light, balanced and the small size mimics a handpiece providing comfortable, clinical familiarity.



Clinically familiar position and technique
Dr. Donald Erickson, DDS, using the 3M™ True Definition Scanner.



The real value of digital technology is in the data and what you can do with it.

1. Study methodology described in: van der Meer WJ, et. al. (2012). Application of Intra-Oral Dental Scanners in the Digital Workflow of Implantology. PLoS ONE 7(8):e43312. doi:10.1371/journal.pone.0043312. Additional measurements conducted by ACTA (Academic Center for Dentistry Amsterdam); Wicher J. van der Meer, et. al. (2012). Publication pending.

3M is boldly redefining the economics of CAD/CAM dentistry

Introducing

3M True Definition Scanner

The 3M™ True Definition Scanner is shattering the high financial barrier of CAD/CAM dentistry with an unprecedented price.

The 3M True Definition system is supported by an open and secure cloud-based platform, offering you and your dental lab more ways to store clinical cases, share and connect.

Finally, there's an advanced digital impression system that's feasible and a smart investment for any dental practice—whether well-established or just starting out.

See why it's time to
go digital with 3M at

www.3M.com/TrueDef



3M ESPE

Dental Products

2510 Conway Avenue
St. Paul, MN 55144-1000 USA
1-800-634-2249

Please recycle. Printed in U.S.A.
70-2013-0536-7

3M is a trademark of 3M or
3M Deutschland GmbH. Used under
license in Canada. © 3M 2012.
All rights reserved.