

QUICKSCANIOS™




The **Full Color, Powderless** Intraoral Scanner

Accurate, Digital Impression Scans

The QuickScan IOS is powderless and offers full color digital scans. It is equipped with the latest technology of the new generation of intraoral scanners. Tested by independent labs, it's accuracy and performance is comparable to other leading scanners but at a fraction of the cost!

COMPATIBLE WITH

 QuickScan IOS scans are powered by Exocad and can be viewed in leading software such as 3shape, Bluesky Bio and more that accept stl and ply files.



Your Practice Deserves To Purchase The *Right* Scanner!
 Let Our Team of Experts Help You Make the Best Choice


TECHNICAL SPECIFICATIONS

SCANNER SIZE	240 x 49.5 x 30.5mm
SCANNER WEIGHT	210g
DISTANCE	0-17mm
ACCURACY	20um
IMAGING MODE	Synchronous 3D Video
QUICK RESET	Yes
MODEL OPTIMIZATION	Yes



FEATURES

- **Lightweight, Ergonomic, Handpiece**
- **Powderless Scanning**
- **Capable of Viewing Scans in Full Color or Monochrome**
- **Full Arch Scanning in Less Than 3 Minutes**
- **Heated Autoclavable Tips to Control Fogging**
- **Can Accomodate all Dental Workflows Including Crowns, Bridges, Implant Planning, Implant Restorations, Complete Dentures, Smile Design & Orthodontics, to Name a Few**
- **Open Architecture: Export Your Scans as STL Files to all Open Treatment Planning & Design Software: Exocad, 3Shape, BlueSkuBio, etc**
- **1 Year Warranty with Optional 3 Year Replacment.**



Produces STL + PLY Files

OPTIONAL ADD ONS



Laptop



Touchscreen Computer Cart



Autoclavable Tips

[DEALER TAG HERE]

CONTACT YOUR DEALER TODAY!

DENTERPRISE
 INTERNATIONAL, INC

Denterpriseintl.com