

STRAUMANN SIRIOS™

Capture Excellence.



STRAUMANN SIRIOS™ **CAPTURE** **EXCELLENCE**


Straumann SIRIOS™ goes beyond image capturing, redefining the essence of how you interact with dentistry and opening the door to excellence, with your patient at the center.

The synergy of our performance-driven hardware with our intuitive software is engineered to deliver peak performance across your workflow with precision and efficiency.






**REVOLUTIONARY
SCANNING PERFORMANCE**



**DESIGNED WITH YOUR
WORKFLOW IN MIND**



**UNLOCK ULTIMATE
ACCURACY**



**ELEVATED STRAUMANN®
EXPERIENCE**

“The Straumann SIRIOS™ gives me the scanning performance I need, while allowing me to quickly incorporate it into my existing workflows. Its wireless design, with optimal size, weight, and handling – makes scanning fast, easy and comfortable for my patients, whilst giving me the flexibility to adapt to individual treatment needs.”

Dr. Gustavo Harfagar

REVOLUTIONARY SCANNING PERFORMANCE

Redefine ease of use with our intraoral scanning solution. Straumann SIRIOS™ features a smart design and enables a flexible scanning strategy, empowering you to adapt to several clinical situations. Start scanning easily to deliver high-speed, optimal results efficiently.

DESIGN

LIGHTWEIGHT AND COMPACT DESIGN

This wireless scanner combines cutting-edge technology with an ergonomic pen-grip design for user and patient comfort. Weighing only 245 g (including battery), it ensures effortless scanning and reduces fatigue during extended use.

TWO TIP SIZES

Select the tip size that best suits the clinical situation to ensure a premium care experience for your patients.

FLEXIBLE SCANNING STRATEGY

Enjoy flexibility in your scanning process. Adapt easily to each clinical situation, without the need for movement restriction or mandatory scanning strategy.





EFFICIENCY

WIRELESS CHARGING

Charging is as easy as placing the device back into its holder, for effortless all-day power. With Instant Wake-up, there is no waiting time, it activates the moment it is removed from the holder.

UP TO TWO HOURS OF CONTINUOUS SCANNING

Enjoy two hours of continuous scanning on a single charge, providing power to handle the equivalent of 60 standard cases, and ensuring you stay efficient throughout your workday.

HIGH-SPEED SCANNING

Complete a full-arch scan in under 20 seconds with high-quality data capture and real-time visualization. Experience ultra-fast, smooth connectivity thanks to our advanced image compression technology.

SIMPLICITY

CALIBRATION FREE

Eliminate the need for costly calibration tools and avoid disrupting your routine. We ensure the components maintain their integrity, making the device ready for precise scanning from day one and reliable over time.

DUO SMART SCANNING BUTTON

Two strategically placed scan buttons on either side of the scanner allow for seamless rotation and continuous scanning without the need to remember button positions or reassemble the tip.

EASY TO INSTALL AND MASTER

Wireless, our solution offers a quick, true plug-and-play setup, allowing you to start scanning immediately. Start scanning effortlessly with our optimized scanning technology and user-friendly software.

DESIGNED WITH YOUR WORKFLOW IN MIND

As the doorway to the Straumann® integrated digital ecosystem, Straumann SIRIOS™ is designed to unlock your workflow to achieve excellence. Optimized for Straumann® solutions, along with integration into the new Straumann AXS™, it provides you with seamless and tailored workflows to respond to your needs.

OPTIMIZED FOR STRAUMANN® SOLUTIONS

Straumann SIRIOS™ is integrated into the Straumann® digital ecosystem with seamless workflow solutions.

IMPLANT WORKFLOW

Automatic scanbody recognition not only streamlines your workflow but also validates the quality of your scan.

Dedicated and validated full-arch workflow to achieve high accuracy for advanced implant restorative treatment, all combined in one single scanning case. Gain access with even superior user experience when using our Straumann EXACT™ solution.

SURGICAL PLANNING WORKFLOW

Plan implant cases and design surgical guides easily. Export IO scans to coDiagnostiX® installed on your laptop.

ORTHODONTIC WORKFLOW

Treat aligner, revision and retainer cases from the Straumann SIRIOS™ software. Have direct patient-case connectivity to Dr. Portal from ClearCorrect®.



STREAMLINED AND PEACE-OF-MIND WORKFLOW POWERED BY STRAUMANN AXS™

Gain access to everything in one place through Straumann AXS™. Your end-to-end digital workflows are now totally seamless.

Experience seamless connectivity and effortless workflow integration with Straumann SIRIOS™ and Straumann AXS™: your IOS scan files are automatically uploaded to Straumann AXS™, ensuring secure storage and easy access to the next steps in treatment. With synchronized patient data between both solutions, avoid duplicate entries and streamline efficiency for optimal treatment outcomes.

SEND TO LAB WORKFLOW:


Streamline your practice and communication, easily sharing cases with your lab in only one click.

SMILE IN A BOX® WORKFLOW:

Outsource guided surgery cases easily. Link IOS scans to the Smile in a Box® case order form in Straumann AXS™.

STRAUMANN UN!Q™ WORKFLOW:

Outsource design and production of patient-specific restorations easily on Straumann AXS™.



Discover more at:
www.straumann.com/axs

UNLOCK ULTIMATE ACCURACY

Empower your clinical decision-making with high predictability, thanks to high-resolution scanning and accurate 3D files. After experiencing our advanced scanning technology and selecting your preferred optimized workflow, unlock numerous treatment possibilities and deliver an enhanced patient experience with Straumann SIRIOS™.

MINIMIZE REWORK

Straumann SIRIOS™ provides precise digital impressions, significantly reducing errors compared to plaster models.

You can rely on scanbody recognition and validation, eliminating errors and streamlining operations and communications between CAD software and lab partners.

SAVE VALUABLE TIME

Straumann SIRIOS™ digital impressions can be stored and easily reused for future restorations. This eliminates the need for new scans saving valuable time and reducing rework costs.

Our AI-driven software helps you save valuable chair time by automating critical scanning steps. Especially in advanced cases, it enhances efficiency by duplicating IOS scans and intelligently filtering out the scanbodies, allowing you to rescan only clinically relevant areas with ease.

OPTIMIZE CLINICAL OUTCOMES

Enables highly accurate scanning, empowering you to achieve optimal clinical outcomes. Enhanced fit and function, combined with the ability to offer a broader range of treatment options, ultimately lead to higher patient satisfaction.

By offering an end-to-end validated solution from the data acquisition to the final prosthesis manufacturing, with Straumann SIRIOS™ and Straumann EXACT™ you can deliver reliable clinical outcomes for full-arch cases.



Software interface subject to updates for continuous improvement.

**Smart Scanning Button**

Start scanning with a single touch and monitor your progress in real-time with dynamic lighting guidance.

**Mode Switch Button**

Easily navigate through different workflow stages with a click.

**Status Indicator**

Check scanning status and progress with lighting indication.

**Battery Indicator**

Be reminded about battery level status with lighting indication.



THE ELEVATED STRAUMANN® EXPERIENCE

From Straumann®, you can expect nothing less than premium service and customer experience that follow you since the beginning of your journey with us.

SUPERIOR TRAINING AND EDUCATION

Our commitment to education is a cornerstone of effective patient care and practice success. We offer a wide range of educational programs, including classroom courses, eLearning, webinars – all integrated into one digital platform – worldwide.

For Straumann SIRIOS™ you can also count on a unique e-Learning Masterclass onboarding. It features dynamic and innovative content to provide you with an engaging learning experience.

SUPPORT YOU CAN RELY ON

Have access to a customer support program that offers quick issue resolution, availability in multiple geographic locations, and local language assistance.

Benefit from a comprehensive customer care package that includes extended warranty, full accidental coverage, and all the valuable Straumann® customer care services you already know and trust.



Straumann® is the flagship brand of the world's largest implant manufacturer group – the Straumann Group. Renowned as one of the most iconic brands in dentistry, we have a long-standing legacy of pioneering significant advancements in the field. We are at the forefront of expanding the scope of dental care by relentlessly and consistently providing excellence and transformational innovation to a global community of forward-thinking dental professionals.



STRAUMANN SIRIOS™ CAPTURE **EXCELLENCE.**



490.950/en/B/00 03/25

International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone +41 (0)61 965 11 11
Fax +41 (0)61 965 11 01
www.straumann.com

© Institut Straumann AG, 2025. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

 **straumann**